



شرکت پارسیان گیتی ایرانیان با هدف تولید و ارائه جدیدترین و با کیفیت ترین محصولات در زمینه تجهیزات پزشکی، آزمایشگاهی در سال ۱۳۸۸ تاسیس گردید. امروزه با تکیه به سابقه و تجربه چهارده ساله شرکت در زمینه واردات و تولید و داشتن گروهی متخصص در کنار تجهیزاتی بسیار پیشرفته در بخش آزمایشگاهی توانسته ایم تبدیل به بزرگترین تولیدکننده در حوزه کیت های تشخیص سریع (رپید تست ها) شویم همچنین در نوآوری، تولید و واردات ملزومات پزشکی، آزمایشگاهی و تجهیزات پزشکی، هم راستا با کمپانی های بزرگ دنیا و استانداردهای جهانی باشیم.

برای معرفی برخی از محصولات شرکت می توان از انواع کیت های تشخیص سریع بارداری و اعتیادی، تست های تخمک گذاری و کیت های تشخیص سریع ویروس کووید و تولید ماسک های سه لایه پزشکی نام برد.

Parsian Giti Iranian Company was established in 2009 with the aim of producing and presenting the latest and highest quality products in the field of medical equipment. Today, relying on the fourteen years of experience and experience of the company in the field of import and production and having a specialized team along with highly advanced equipment in the laboratory, we have been able to become the largest manufacturer in the field of rapid detection kits (rapid tests) Import of medical supplies, laboratory and medical equipment, in line with the world's major companies and international standards.

To introduce some of the company's products, we can mention various types of rapid pregnancy and addiction diagnostic kits, ovulation tests and rapid detection kits of Covid virus and production of three-layer medical masks.



PREGNANCY TEST

Medicore Test
more Care. more Health

MEDI  **HEALTH**
science & quality

A pregnancy test can tell whether you are pregnant by checking for a particular hormone in your urine or blood. The hormone is called human chorionic gonadotropin (HCG). HCG is made in a woman's placenta after a fertilized egg implants in the uterus. It is normally made only during pregnancy



Pregnancy Test | Cassette

- Cut off : 25 ng/mL
- Accuracy : 99.8%
the first morning urine usually contains the highest amount of hCG that is preferred.
- Detection of pregnancy before expected period.
- The most widely used pregnancy diagnostic test in country
- Market value, Convenient and easy to use, With dropper
- Licensed by the department of medical equipment



Quality Management System
EN ISO 13485 : 2016 / Made in IRAN

 **Parsian giti iranian Co., Ltd**
Km10 Karaj-Ghazvin Hw. Barkat Zone

www.parsiangiti.co



Early Pregnancy Detection Matters

The first few weeks of pregnancy are critical to fetal development, so early prenatal care is important for a healthy birth weight and the baby's survival. The sooner you know you're pregnant, the sooner you can begin making healthy pregnancy choices for yourself and your baby that include diet, nutrition and limiting your exposure to hazards.



Pregnancy Test | Strip

- Cut off : 25 ng/mL
- Accuracy : 99.8%
the first morning urine usually contains the highest amount of hCG that is preferred.
- Detection of pregnancy before expected period.
- The Most widely used diagnostic test in the country
- Market value
- Convenient and easy to use
- Licensed by the department of medical equipment



Quality Management System
EN ISO 13485 : 2016 / Made in IRAN



Parsian giti iranian Co.,Ltd
Km10 Karaj-Ghazvin Hw. Barkat Zone

www.parsiangiti.co



MEDI HEALTH[®]
science & quality Test

تست بارداری | نواری

• Cut off : 25 ng/mL

• بالاترین دقت : ۹۹.۸ درصد

بهترین نمونه ادرار ، اولین ادرار صبح به علت غلظت بالای آن است

• تشخیص بارداری قبل از شروع دوره قاعدگی

• پر مصرف ترین تست تشخیص بارداری در کشور

• قیمت مناسب، کاربرد بسیار راحت

• دارای مجوز از اداره کل تجهیزات و ملزومات پزشکی

Quality Management System
EN ISO 13485 : 2016 / Made in IRAN



Parsian giti iranian Co.,Ltd
Km10 Karaj-Ghazvin Hw. Barkat Zone

www.parsiangiti.co



Pregnancy Test | Midstream

- Cut off : 25 ng/mL
- Accuracy : 99.8%
the first morning urine usually contains the highest amount of hCG that is preferred.
- Detection of pregnancy before expected period.
- Best use in vacations
- Market value, Convenient and easy to use
- Licensed by the department of medical equipment



Quality Management System
EN ISO 13485 : 2016 / Made in IRAN



Parsian giti iranian Co., Ltd
Km10 Karaj-Ghazvin Hw. Barkat Zone

www.parsiangiti.co



تست بارداری | خودکاری

Cut off : 25 ng/mL •

• بالاترین دقت : ۹۹.۸ درصد

بهترین نمونه ادرار ، اولین ادرار صبح به علت غلظت بالای آن است

• تشخیص بارداری قبل از شروع دوره قاعدگی

• استفاده راحت در سفر

• قیمت مناسب، کاربرد بسیار راحت

• دارای مجوز از اداره کل تجهیزات و ملزومات پزشکی





Pregnancy Test | Strip

- Cut off : 25 ng/mL
- Accuracy : 99.8%
the first morning urine usually contains the highest amount of hCG that is preferred.
- Detection of pregnancy before expected period.
- Sensitivity = 40 mIU
- Best use in vacations
- Market value, Convenient and easy to use
- Licensed by the department of medical equipment





Medicore[®] Test

more Care , more Health

تست تخمک گذاری | نواری

Cut off : 25 ng/mL •

• بالاترین دقت : ۹۹.۸ درصد

بهترین زمان انجام تست ادرار میانه روز ۱۰ صبح تا ۴ بعداز ظهر می باشد و ترجیحا از ادرار صبحگاهی استفاده نشود

• اطلاع از بهترین زمان برای بارداری

• دقیق ترین تست تشخیص زمان تخمک گذاری در کشور

• کاربرد راحت، دقت بالا

• دارای مجوز از اداره کل تجهیزات و ملزومات پزشکی

Quality Management System
EN ISO 13485 : 2016 / Made in IRAN



Parsian giti Iranian Co.,Ltd
Km10 Karaj-Ghazvin Hw. Barkat Zone

www.parsiangiti.co

SARS-COV 2 ANTIGEN RAPID TEST CASSETTE

FIRST
CHECK[®]

Diabase[®]

Medicore[®]
more Care, more Health

MEDI[®] HEALTH[®]
science & quality

Antigen tests are immunoassays that detect the presence of a specific viral antigen, which implies current viral infection. Antigen tests are currently authorized to be performed on nasopharyngeal or nasal swab specimens placed directly into the assay's extraction buffer or reagent. The currently authorized antigen tests are not restricted to use on persons of a certain age.

Diabase[®] ±

SARS - COV-2 Antigen | Rapid Test Cassette

- Cut off : 50 ng/mL
- Sensitivity: 97.9%
- Screening through a sample of Nasopharynx swap
- Answer in 15 minutes
- Pretty convenient to use, reasonable price
- Has nasopharyngeal swab & buffer extraction solution
- Licensed by the department of medical equipment

15 Just
Minutes

Highest
Quality

Quality Management System
EN ISO 13485 : 2016 / Made in IRAN



Parsian giti iranian Co., Ltd
Km10 Karaj-Ghazvin Hw. Barkat Zone

www.parsiangiti.co





Diabase[®] ±

تست سریع آنتی ژن سارس - کووید ۲ | کاسته

- Cut off : 50 ng/mL
- غربالگری از طریق نمونه سواب نازوفارنکس
- حساسیت % 97.9 |
- جواب دهی در ۱۵ دقیقه
- کاربرد بسیار راحت ، قیمت مناسب
- دارای سواب نازوفارنکس و محلول استخراج بافر
- دارای مجوز از اداره کل تجهیزات و ملزومات پزشکی

Quality Management System
EN ISO 13485 : 2016 / Made in IRAN



Parsian giti Iranian Co.,Ltd
Km10 Karaj-Ghazvin Hw. Barkat Zone

www.parsiangiti.co



SARS - COV-2 Antigen | Rapid Test Cassette

- Cut off : 50 ng/mL
- Sensitivity: 97.9%
- Screening through a sample of Nasopharynx swap
- Answer in 15 minutes
- Pretty convenient to use, reasonable price
- Has nasopharyngeal swab & buffer extraction solution
- Licensed by the department of medical equipment



Quality Management System
EN ISO 13485 : 2016 / Made in IRAN

 Parsian giti iran Co.,Ltd
Km10 Karaj-Ghazvin Hw. Barkat Zone

www.parsiangiti.co





FIRST[®] CHECK

تست سریع آنتیژن سارس - کووید ۲ | کاستی

- Cut off : 50 ng/mL
- غربالگری از طریق نمونه سواب نازوفارنکس
- حساسیت %97.9
- جواب دهی در ۱۵ دقیقه
- کاربرد بسیار راحت ، قیمت مناسب
- دارای سواب نازوفارنکس و محلول استخراج بافر
- دارای مجوز از اداره کل تجهیزات و ملزومات پزشکی



Quality Management System
EN ISO 13485 : 2016 / Made in IRAN




Parsian giti Iranian Co.,Ltd
Km10 Karaj-Ghazvin Hw. Barkat Zone

www.parsiangiti.co

SARS - COV-2 Antigen | Rapid Test Cassette

- Cut off : 50 ng/mL
- Sensitivity: 97.9%
- Screening through a sample of Nasopharynx swap
- Answer in 15 minutes
- Pretty convenient to use, reasonable price
- Has nasopharyngeal swab & buffer extraction solution
- Licensed by the department of medical equipment





Medicore[®] Test

more Care , more Health

تست سریع آنتیژن سارس - کووید ۲ | کاستی

- Cut off : 50 ng/mL
- غربالگری از طریق نمونه سوآپ نازوفارنکس
- حساسیت %97.9
- جواب دهی در ۱۵ دقیقه
- کاربرد بسیار راحت ، قیمت مناسب
- دارای سوآپ نازوفارنکس و محلول استخراج بافر
- دارای مجوز از اداره کل تجهیزات و ملزومات پزشکی



Quality Management System
EN ISO 13485 : 2016 / Made in IRAN



Parsian giti Iranian Co.,Ltd
Km10 Karaj-Ghazvin Hw. Barkat Zone

www.parsiangiti.co

SURGICAL FACE MASK

Medicore[®]
Test
more Care . more Health

MEDI  **HEALTH**[®]
science & quality

Social distancing is one of the most important things you can do during a pandemic to help save lives and prevent our health care system from exceeding capacity. So, when you do go out of the house, be sure you continue to practice social distancing and wear a cloth face cover when going into businesses or into gatherings other than with members of your household

Surgical Face Mask



Mask | 3ply surgical face mask

Ear loop & Tie-on high filtering mask

- Has a Melt Blown filter with medical grade
- Filtration of particles to 0.25 micron
- Hypoallergenic
- Very low resistance to breathing
- It has a ergonomic structure and is easy to use
- Licensed by the department of medical equipment
- Approved by Nuclear Science and Technology Research Institute



*Fresh
Breathable*

Quality Management System
EN ISO 13485:2016 / Made in IRAN

 **Parsian giti iran Co., Ltd**
Km10 Karaj-Ghazvin Hw. Barkat Zone

www.parsiangiti.co





MEDI HEALTH
science & quality [®] Test

ماسک | جراحی ۳ لایه ماسک فیلتردار کش دار و بنددار

- حاوی فیلتر Melt Blown با گرید پزشکی
- قابلیت فیلتراسیون ذرات معلق تا ۰.۲۵ میکرون
- مقاومت بسیار پایین در مقابل تنفس، ضد آلرژی
- دارای ساختار ارگونومیک و کاربرد آسان
- دارای مجوز از اداره کل تجهیزات و ملزومات پزشکی
- دارای تاییدیه از پژوهشگاه علوم و فنون هسته ای سازمان انرژی اتمی

Masks are recommended as a simple barrier to help prevent respiratory droplets from traveling into The air and onto other people when the person wearing the mask coughs, sneezes, talks, or raises their voice.

Face Masks We Actually Like to Wear We'll be wearing them for a while, so here are the favorite cloth coverings for running, walking, and going to work

Quality Management System
EN ISO 13485 : 2016 / Made in IRAN

Parsian giti iran Co.,Ltd
Km10 Karaj-Ghazvin Hw. Barkat Zone

www.parsiangiti.co

Mask | 3ply surgical face mask

Ear loop & Tie-on high filtering mask

- Has a Melt Blown filter with medical grade
- Filtration of particles to 0.25 micron
- Hypoallergenic
- Very low resistance to breathing
- It has a ergonomic structure and is easy to use
- Licensed by the department of medical equipment
- Approved by Nuclear Science and Technology Research Institute

FAMILY PROTECTED

Social distancing is one of the most important things you can do during a pandemic to help save lives and prevent our health care system from exceeding capacity. So, when you do go out of the house, be sure you continue to practice social distancing and wear a cloth face cover when going into businesses or into gatherings other than with members of your household

ماسک | جراحی ۳ لایه

ماسک فیلتردار کش دار و بنددار

- حاوی فیلتر Melt Blown با گرید پزشکی
- قابلیت فیلتراسیون ذرات معلق تا ۰/۲۵ میکرون
- پنج عدد ماسک داخل یک کاور سلفونی
- مقاومت بسیار پایین در مقابل تنفس ، ضد آلرژی
- دارای ساختار ارگونومیک و کاربرد آسان
- دارای مجوز از اداره کل تجهیزات و ملزومات پزشکی
- دارای تاییدیه از پژوهشگاه علوم و فنون هسته ای سازمان انرژی اتمی



We Protect You!

Mask | 3ply surgical face mask

Masks are recommended as a simple barrier to help prevent respiratory droplets from traveling into the air and onto other people when the person wearing the mask coughs, sneezes, talks, or raises their voice. This recommendation is based on what we know about the role respiratory droplets play in the spread of the virus that causes COVID19-, paired with emerging evidence from clinical and laboratory studies that shows masks reduce the spray of droplets when worn over the nose and mouth.



Quality Management System
EN ISO 13485 : 2016 / Made in IRAN



Parsiangiti Iranian Co., Ltd
Km10 Karaj-Ghazvin Hwy. Barkat Zone

www.parsiangiti.co



DRUG OF ABUSE

Medicore[®]
more Care. more Health

MEDI  **HEALTH**[®]
science & quality

People from all walks of life can experience problems with their drug use, regardless of age, race, background, or the reason they started using drugs in the first place. Some people experiment with recreational drugs out of curiosity, to have a good time, because friends are doing it, or to ease problems such as stress, anxiety, or depression.



Drug Abuse | AMPHETAMINE, METHAMPHETAMINE

- Cut off : 500 ng/mL
- Accuracy : 99.9%
- Urine samples can be taken at any time of the day, Answer possibility of up to 72 hours after consumption
- Test Used for each group of glass or crystal
- Convenient and easy to use, answers in 3 to 5 minutes



تست‌های اعتیادی | آمفتامین مت‌آمفتامین

- Cut off : 300 ng/mL
- بالاترین دقت : ۹۹.۹ درصد
- امکان نمونه گیری ادرار در هر ساعت از روز
- جواب دهی تا ۷۲ ساعت بعد از مصرف
- قابل استفاده برای هر گروه شیشه یا کریستال
- کاربرد راحت، جواب دهی ۳ الی ۵ دقیقه



100 Tests



Quality Management System
EN ISO 13485 : 2016 / Made in IRAN

 **Parsian giti Iranian Co., Ltd**
Km10 Karaj-Ghazvin Hw. Barkat Zone

www.parsiangiti.co

Drug Abuse | MORPHINE , METHAMPHETAMINE , MARIJUANA - CANNABIS

- Cut off MOP : 300 ng/mL , Cut off MET : 500 ng/mL , Cut off THC : 50ng/ml
- Accuracy : 99.9%
- Urine samples can be taken at any time of the day, Answer
- possibility of up to 72 hours after consumption
- Reducing time & cost
- Convenient and easy to use, answers in 3 to 5 minutes
- Licensed by the department of medical equipment

تست‌های اعتیادی | مورفین ، متادون ، حشیش

- بالاترین دقت : ۹۹.۹ درصد
- امکان نمونه گیری ادرار در هر ساعت از روز
- جواب دهی تا ۷۲ ساعت بعد از مصرف
- کاهش هزینه و زمان
- کاربرد راحت، جواب دهی ۳ الی ۵ دقیقه
- دارای مجوز از اداره کل تجهیزات و ملزومات پزشکی
- بالاترین دقت : ۹۹.۹ درصد



100 Tests

Quality Management System
EN ISO 13485 : 2016 / Made in IRAN



Parsian giti iranien Co.,Ltd
Km10 Karaj-Ghazvin Hw. Barkat Zone

www.parsiangiti.co